Oxy₂Mask™

O₂ / EtCO₂

REF

Oxygen Tube Length	7 feet		14 feet	
Sampling Line	SLF 8"	SLM 8'	SLF 8"	SLM 15'
Oxy₂Mask™ EtCO₂ Kid	OK-2125-8	OK-2125-8SLM	-	-
Oxy₂Mask™ EtCO₂ Adult	OM-2125-8	OM-2125-8SLM	OM-2125-14	OM-2125-14SLM
Oxy₂Mask™ EtCO₂ Plus	OP-2125-8	OP-2125-8SLM	OP-2125-14	OP-2125-14SLM

Intended Use: A minimal contact, open oxygen delivery device for use with USP medical grade oxygen for monitoring breathing by providing a means to sample exhaled CO₂. **Warning:** Hypercapnia will not be reliably detected by sampling through this device, the displayed CO₂ values or any other parameter/index driven from CO₂ will vary with the sampling port position, oxygen flow rate and respiratory pattern.

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Caution: Oxygen in the presence of an ignitable material and an ignition source can cause a fire resulting in severe personal injury or death. Do not smoke while using this device. Single use product not to be reused due to effect on performance, cross contamination and risk of infection. The device has a 5-year shelf life from the date of manufacture. Date of manufacture is stated on labelling. Store device in original packaging until use. Inspect product and packaging prior to use. Do not use if product appears to be damaged.

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
EC REP
EMERGO EUROPE
Westervoortsejdijk 60
6827 AT Arnhem
The Netherlands

CH REP
MedEnvoy Switzerland
Gotthardstrasse 28
6302 Zug
Switzerland

UK CA 0086 **CE** 2797
Emergo Consulting (UK) Limited
c/o Cr 360—UL International
Compass House, Vision Park Histon
Cambridge, CB24 9BZ
United Kingdom

IFU0299E Issue Date: 2024-03-06

50 Alliance Blvd. Barrie, Ontario Canada L4M 5K3 705-726-9383 705-728-9537 FAX thebetteroxygenmask.com

EC REP
EMERGO EUROPE
Westervoortsejdijk 60
6827 AT Arnhem
The Netherlands

CH REP
MedEnvoy Switzerland
Gotthardstrasse 28
6302 Zug
Switzerland

UK CA 0086 **CE** 2797
Emergo Consulting (UK) Limited
c/o Cr 360—UL International
Compass House, Vision Park Histon
Cambridge, CB24 9BZ
United Kingdom

IFU0299E Issue Date: 2024-03-06